Ceramics

Analysis of kaolinite

Experimental:

Instrument: Setsys TGA with DTA rod Pt-Rh 10% / Pt Sample: kaolinite Sample mass: 78.87 mg Crucible: platinum Atmosphere: air The temperature is programmed from 25°C up to 1600°C at 10 K.min⁻¹. Notice: kaolinite is the raw material for the production of china (porcelain).



Conclusion:

The TG curve presents two mass losses:

- the first mass loss of 4.81% between the ambiant and 250° , simultaneous with an endotherm, is due to the evaporation of water.

- the second mass loss of 12.22% between 250°C and 940°C, simultaneous with an endotherm, is due to the dehydroxylation of kaolinite:

kaolinite \rightarrow metakaolinite

The exotherm at 993℃ is to the transformation:

metakaolinite \rightarrow mullite

The exotherm at 1263° is to the transformation: mullite \rightarrow secondary mullite and cristobalite

> Instrument: Setsys Evolution TGA-DTA -150℃ to 2400℃





Excellence in thermal analysis and calorimetry